

Super Sensitive IHC Detection System Kit (Mouse/Rabbit)

Cat No: BD6001

introduction:

Super Sensitive IHC Detection System Kit is the latest technology in polymeric labeling. Polymer detection methods have been shown to provide increased sensitivity. This innovative polymer technology has major advantages than conventional IHC systems. Super Sensitive IHC Detection System amplifies the signal with both mouse and rabbit primary antibodies. Background noise due to nonspecific binding to endogenous biotin molecules is eliminated, because Super Sensitive IHC Detection System is not a biotin/avidin based system, which eliminates the background noise due to nonspecific binding to endogenous biotin. The Super Sensitive IHC Detection System Kit provides the user with a rapid, easy to use, and versatile IHC detection system.

Reagents:

- A: Hydrogen Peroxide Blocking Reagent 5 ml
- B: Blocking Reagent 5 ml
- C: HRP Polymer 5 ml
- D: DAB Substrate 5 ml
- E: DAB Chromogen 100 μ l
- F: Technical Manual 1 Manual

Storage&Shelf life:

All kit components are stable at 4 °C. Each component is stable for up to 12 months.

Procedure:

1. Deparaffinize and rehydrate tissue section; PBS/TBS wash 3 times for 2 minutes;
2. Incubate tissue in appropriate pretreatment or digestive enzyme if required for primary antibody; and PBS/TBS wash 3 times for 2 minutes;
3. Incubate slide in Hydrogen Peroxide Blocking Reagent for 10 minutes, PBS/TBS wash 3 times for 2 minutes;
4. Apply Blocking Reagent and incubate for 5 minutes, PBS/TBS wash 3 times for 2 minutes (May be omitted if


primary antibodies are diluted in buffers containing normal goat serum);

5. Apply HRP Polymer (50 μ l for each slice) and incubate for 30 minutes, PBS/TBS wash 3 times for 2 minutes;

6. Add 20 μ l DAB Chromogen into 1 ml DAB Substrate, mixing vial shortly before incubation, and apply to tissue (50 μ l for each slice). Incubate for about 3 - 5 minutes, PBS/TBS wash for 2 minutes;

7. Counterstain and coverslip using a permanent mounting media.

DATA:



Immunohistochemical analysis staining in human breast carcinoma formalin fixed paraffin-embedded tissue section. The section was pre-treated using pressure cooker heat antigen retrieval with sodium citrate buffer (0.01M, pH=6) for 3 minutes. The section was detected using Super Sensitive IHC Detection System. The section was then counterstained with haematoxylin and mounted with Neutral Gum.

Research Use:

For research use only. Not for diagnostic or therapeutic procedures.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416, USA.

Email: info@bioworld.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046, P. R. China.

Email: info@biogot.com

Tel: 0086-025-68037686

Fax: 0086-025-68035151