

PRODUCT DATA SHEET

Bioworld Technology, Inc.

AMPKα1 polyclonal antibody

Catalog: BS90052 Host: Rabbit Reactivity: Human, Mouse, Rat

BackGround:

AMPK (for 5'-AMP-activated protein kinase) is a heterotrimeric complex com-prising a catalytic α subunit and regulatory β and γ subunits. It protects cells from stresses ATP depletion by switching ATP-consuming bio-synthetic pathways. AMPK is activated by high AMP and low ATP through a mechanism involving allosteric regulation, promotion of phosphorylation by an upstream protein kinase known as AMPK kinase, and inhibition of dephosphorylation. Activated AMPK can phosphorylate and regulate in vivo hydroxy-methylglutaryl-CoA reductase and acetyl-CoA carboxylase, which are key regulatory enzymes of sterol synthesis and fatty acid synthesis, respectively. The human AMPKα1 and AMPKα2 genes encode 548 amino acid and 552 amino acid proteins, respectively. Human AMPKβ1 encodes a 271 amino acid protein and human AMPKβ2 encodes a 272 amino acid protein. The human AMPKy1 gene encodes a 331 amino acid protein. Human AMPKγ2 and AMPKγ3, which are 569 and 492 amino acid proteins, respectively, contain unique N-terminal domains and may participate directly in the binding of AMP within the AMPK complex.

Product:

Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

Molecular Weight:

63 kDa

Swiss-Prot:

Q13131(Human) Q5EG47(Mouse) P54645(Rat)

Purification&Purity:

ProA affinity purified

Applications:

WB:1:1,000-1:5,000 ICC:1:50-1:200

FC:1:50-1:100


Storage&Stability:

Store at +4 $^{\circ}$ C after thawing. Aliquot store at -20 $^{\circ}$ C or -80 $^{\circ}$ C. Avoid repeated freeze / thaw cycles.


Specificity:

AMPK α 1 polyclonal antibody detects endogenous levels of AMPK α 1 protein.

DATA:



Western blot analysis of AMPK alpha 1 on different lysates using anti-AMPK alpha 1 antibody at 1/1,000 dilution. Positive control: Lane 1: Hela Lane 2: HepG2 Lane 3: MCF-7



ICC staining AMPK alpha 1 in Hela cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: <u>info@bioworlde.com</u>

Tel: 6123263284 Fax: 6122933841

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: <u>info@biogot.com</u>
Tel: 0086-025-68037686
Fax: 0086-025-68035151