

PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



WB Blocking Buffer

Cat No.: WBA2400-R1

Introduction

This product is used to block PVDF or NC membranes. The main active can effectively block the protein binding sites on membranes and reduce non-specific bindings of antibodies, thus reducing the background and enhancing the signal-to-noise ratio.

Quantity

100ml/500mL

Procedure

This product is a ready-to-use liquid and can be directly used for blocking Western Blot membranes. Usually, blocking for 60 minutes at room temperature is sufficient.

Storage & Shelf life

Store at room temperature. Each component is stable for up to 6 Months.

NOTE

If the buffer color and smell is changed, you'd better to discard it.

Research Use

For research use only, not for use in diagnostic procedures.

MADE IN CHINA

Bioworld Technology, Inc.
1660 South Highway 100, Suite 500 St. Louis Park, MN
55416, USA. Email: info@bioworld.com
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.
No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com
Tel: 0086-025-68037686 Fax: 0086-025-68035151