

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



# Eosin Y Staining Buffer

Cat No.: BD5036

### Introduction

Eosin Y is a stain with histological applications. This formula is based on the original formulation developed by our laboratory and requires longer application periods. Constituents: Eosin Y and Alcohol. The cytoplasm is stained with pink color.

### Reagents

Eosin Y Staining Buffer 500ml

### Application

For cytoplasm staining.

### Storage & Shelf life

Store at RT. Each component is stable for up to 12 months.

### Procedure

1. Rehydrate tissue if needed;
2. Stain cytoplasm 10 seconds minutes with Eosin (the time depended on section thickness);
3. Rinse in 95% alcohol for 13 seconds (the time depended on section thickness);
4. Other experiments

### Research Use

For research use only, not for use in diagnostic procedures.

**MADE IN CHINA**

**Bioworld Technology, Inc.**

**1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA.**

**Email: [info@bioworld.com](mailto:info@bioworld.com)**

**Tel: 6123263284**

**Fax: 6122933841**

**Bioworld technology, co, Ltd.**

**No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China.**

**Email: [info@biogot.com](mailto:info@biogot.com)**

**Tel: 0086-025-68037686**

**Fax: 0086-025-68035151**