

PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



Neutral Gum

Cat No.: BD5044

Introduction

This product is used for mounting in IHC test. When it becomes density, you can dilute it with xylol.

Reagents

Neutral Gum 100ml

Application

Used for Immunohistochemistry (IHC) experiments.

NOTE

1. Neutral Gum easily volatilizes organic solvents;
2. Operate it in the fuming cupboard;
3. Keep away from fire

Research Use

For research use only, not for use in diagnostic procedures.

MADE IN CHINA

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN
55416, USA. Email: info@bioworlde.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com

Tel: 0086-025-68037686 Fax:0086-025-68035151