

# PRODUCT DATA SHEET

Bioworld Technology, Inc.

# **DSPP** polyclonal antibody

Catalog: BS80527 Host: Rabbit Reactivity: Human, Mouse, Rat

### **BackGround:**

This gene encodes a member of the small integrin-binding ligand N-linked glycoprotein (SIBLING) family of proteins. The encoded preproprotein is secreted by odontoblasts and proteolytically processed to generate two principal proteins of the dentin extracellular matrix of the tooth, dentin sialoprotein and dentin phosphoprotein. These two protein products may play distinct but related roles in dentin mineralization. Mutations in this gene are associated with dentinogenesis imperfecta and dentin dysplasia. This gene is present in a gene cluster on chromosome 4. Allelic differences due to repeat polymorphisms have been found for this gene.

#### **Product:**

1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

### **Molecular Weight:**

110kDa

## **Swiss-Prot:**

Q9NZW4

## **Purification&Purity:**

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

## **Applications:**

WB,1:500 - 1:2000|IF/ICC,1:50 - 1:200

# Storage&Stability:

Store at  $4 \, \mathbb{C}$  short term. Aliquot and store at  $-20 \, \mathbb{C}$  long term. Avoid freeze-thaw cycles.

#### **Modification:**

Unmodification

#### **DATA:**

Western blot analysis of extracts of various cell lines, using DSPP antibody at 1:1000 dilution.<br/>
Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution.<br/>
Lysates/proteins: 25ug per lane.<br/>
Blocking buffer: 3% nonfat dry milk in TBST.<br/>
Detection: ECL Basic Kit.<br/>
Exposure time: 30s.
Immunofluorescence analysis of H9C2 cells using DSPP Rabbit pAb at dilution of 1:100 . Blue: DAPI for nuclear staining.
Immunofluorescence analysis of L929 cells using DSPP Rabbit pAb at

dilution of 1:100 . Blue: DAPI for nuclear staining.

Immunofluorescence analysis of U-2 OS cells using DSPP Rabbit pAb

at dilution of 1:100. Blue: DAPI for nuclear staining.

## Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: <u>info@bioworlde.com</u>

Tel: 6123263284 Fax: 6122933841 Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: <u>info@biogot.com</u>
Tel: 0086-025-68037686
Fax: 0086-025-68035151