

PRODUCT DATA SHEET

Bioworld Technology,Inc.

Thyroid Peroxidase monoclonal antibody

Catalog: MB66022 Host: Mouse Reactivity: Human

BackGround:

The synthesis of thyroid hormones is an oxidative process that produces reac tive oxygen species and requires Thyroperoxidase (TPO), a hemoprotein that is one of the major autoantigens involved in autoimmune thyroid diseases. Thyroperoxidase is a 933 amino acid, type I transmembrane glycoprotein that plays a key role in thyroid hormone synthesis and autoimmunity. TPO catalyzes the iodination of proteins, therefore causing iodide retention within thyroid cells. The ecto-domain of Thyroperoxidase includes a large N-termi nal myeloperoxidase-like domain, followed by a complement control protein domain and an epidermal growth factor-like domain. Thyroperoxidase also mediates the organification and intracellular retention of radioiodide, which may lead to rapid tumor cell death. Mutations of the Thyroperoxidase gene commonly lead to goitrous congenital hypothyroidism, the most severe and frequent abnormality in thyroid iodide organification defect (IOD), in which iodide in the thyroid gland cannot be oxidized and/or bound to the protein.

Product:

Mouse IgG2b. Liquid in PBS containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Molecular Weight:

~ 102 kDa

Swiss-Prot:

P07202

Purification&Purity:

The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

Applications:

IHC (1/100 - 1/300)

Storage&Stability:


Store at 4 ℃ short term. Aliquot and store at -20 ℃ long

term. Avoid freeze-thaw cycles.


Specificity:

Recognizes endogenous levels of Thyroid Peroxidase protein.

DATA:



Immunohistochemical analysis of Thyroid Peroxidase staining in human follicular thyroid carcinoma formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Immunohistochemical analysis of Thyroid Peroxidase staining in human papillary thyroid carcinoma formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: info@bioworlde.com

Tel: 6123263284 Fax: 6122933841

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: <u>info@biogot.com</u> Tel: 0086-025-68037686 Fax: 0086-025-68035151



PRODUCT DATA SHEET

Bioworld Technology,Inc.

Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: <u>info@bioworlde.com</u>
Tel: 6123263284

Fax: 6123263284

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: info@biogot.com
Tel: 0086-025-68037686
Fax: 0086-025-68035151