

# eNOS monoclonal antibody

Catalog: MB65920

Host: N

Mouse

Reactivity: Human, Mouse, Rat

# **BackGround:**

Endothelial nitric-oxide synthase (eNOS) is an important enzyme in the cardiovascular system. It catalyzes the production of nitric oxide (NO), a key regulator of blood pressure, vascular remodeling, and angiogenesis. The activity of eNOS is regulated by phosphorylation at multiple sites. The two most thoroughly studied sites are the activation site Ser1177 and the inhibitory site Thr495. Several protein kinases including Akt/PKB, PKA, and AMPK activate eNOS by phosphorylating Ser1177 in response to various stimuli. In contrast, bradykinin and H2O2 activate eNOS activity by promoting both Ser1177 phosphorylation and Thr495 dephosphorylation.

## **Product:**

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

**Molecular Weight:** 

~ 133 kDa

**Swiss-Prot:** 

# P29474

**Purification&Purity:** 

The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).


## **Applications:**

WB (1/500 - 1/2000), IHC (1/200 - 1/500)


# Storage&Stability:

Store at 4  ${}^\circ\!\!{\rm C}$  short term. Aliquot and store at -20  ${}^\circ\!\!{\rm C}$  long term. Avoid freeze-thaw cycles.


# **Specificity:**



#### **DATA:**



Western blot analysis of eNOS expression in rat heart (A) whole cell lysates.



Immunohistochemical analysis of eNOS staining in mouse colon formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

#### Note:

For research use only, not for use in diagnostic procedure.

# Bioworld Technology, Inc.

 Add:
 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416,USA.

 Email:
 info@bioworlde.com

 Tel:
 6123263284

 Fax:
 6122933841

#### **Bioworld technology, co. Ltd.** Add: No 9, weidi road Qixia District Nanjing, 210046,

 Add:
 Add:
 Add:
 For the formation of the formation o