

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



# Goat Serum Working Buffer

Cat No.: BD5015

### Introduction

This is working goat serum on a certain density. It contains goat serum, PBS, stabilizer and preservative, without Triton X-100. It was aseptically filtered through a Millipore 0.22 um filter into the pre-sterilized container. It is applied for blocking, Immunohistochemistry (IHC), Immunocytochemistry (ICC) and Immunofluorescence (IF).

### Reagents

Goat Serum Working Buffer 100ml

### Application

For blocking in IHC, ICC and IF.

### Storage & Shelf life

Store at 2-8° C for short time; -20°C for long time. Each component is stable for up to 12 Months.

### Research Use

For research use only, not for use in diagnostic procedures.

**MADE IN CHINA**

Bioworld Technology, Inc.  
1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA. Email: [info@bioworld.com](mailto:info@bioworld.com)  
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.  
No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China. Email: [info@biogot.com](mailto:info@biogot.com)  
Tel: 0086-025-68037686 Fax: 0086-025-68035151