## PRODUCT DATA SHEET



**3% Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>)** 

Cat No.: BD5024

#### Introduction

This product is 3% (V/V) hydrogen peroxide. It has been filtered through a 0.22um membrane. It is suitable for e ndogenous peroxydase blocking, DAB and AEC color rea ction.

#### Reagents

3% Hydrogen Peroxide 100ml

#### Application

Blocking for Immunohistochemistry; Substrate Reagent for DAB.

#### **Storage & Shelf life**

Store at  $2-8^{\circ}$  C. Each component is stable for up to 12 Months.

#### NOTE

1. H2O2 has causticity, please operate with gloves and working clothes; keep away from skin;

2. H2O2 resolves out O2, keep away from fire.

#### **Research Use**

For research use only, not for use in diagnostic procedures.

Bioworld Technology, Inc. 1660 South Highway 100, Suite 500 St. Louis Park,MN 55416,USA. Email: <u>info@bioworlde.com</u> Tel: 6123263284 Fax: 6122933841 Bioworld technology, co, Ltd. No 9, weidi road Qixia District Nanjing, 210046, P, R.China. Email: <u>info@biogot.com</u> Tel: 0086-025-68037686 Fax:0086-025-68035151

### MADE IN CHINA

# Bioworld Technology CO., Ltd.