

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



# 30% Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>)

Cat No.: BD5025

### Introduction

This product is 30% (V/V) hydrogen peroxide. It has been filtered through a 0.22um membrane. It is suitable for endogenous peroxidase blocking, DAB and AEC color reaction.

### Reagents

30% Hydrogen Peroxide      30ml

### Application

Blocking for Immunohistochemistry; Substrate Reagent for DAB.

### Storage & Shelf life

Store at 2-8° C. Each component is stable for up to 12 Months.

### NOTE

1. H<sub>2</sub>O<sub>2</sub> has causticity, please operate with gloves and working clothes; keep away from skin;
2. H<sub>2</sub>O<sub>2</sub> releases out O<sub>2</sub>, keep away from fire.

### Research Use

For research use only, not for use in diagnostic procedures.

**MADE IN CHINA**

Bioworld Technology, Inc.  
1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA.      Email: [info@bioworld.com](mailto:info@bioworld.com)  
Tel: 6123263284      Fax: 6122933841

Bioworld technology, co, Ltd.  
No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China.      Email: [info@biogot.com](mailto:info@biogot.com)  
Tel: 0086-025-68037686      Fax: 0086-025-68035151