# **Bioworld Technology CO., Ltd.**



Cat No.: BD5026

# Introduction

This product is 20\*Phosphate Buffered Saline chromogen (0.2M); it contains NaCl, Na<sub>2</sub>HPO<sub>4</sub> and NaH<sub>2</sub>PO<sub>4</sub>. The pH=7.4 $\pm$ 0.1 (0.01M), when it is dissolved in purify water (chromogen: purify water =1:19). This product is a germ free and filtered throughout 0.22um membrane. The dilution buffer is used for Immunohistochemisty (IHC), Immunofluorescence (IF) and so on. It has not been detected for cell culture.

### Reagents

PBS Chromogen 250ml

Application

For washing in IHC, ICC, ELISA, WB and IF.

**Storage & Shelf life** 

Store at  $2-8^{\circ}$  C. Each component is stable for up to 12 Months.

## **Procedure**

**PBS Chromogen** 

Dissolve it in purify water (chromogen: purify water =1:19) for use.

# NOTE

There may some crystalline solid in the buffer when it stored at a low temperature. Take it at RT to make sure there is no crystalline solid when use it.

### **Research Use**

For research use only, not for use in diagnostic procedures.

Bioworld Technology, Inc. 1660 South Highway 100, Suite 500 St. Louis Park,MN 55416,USA. Email: <u>info@bioworlde.com</u> Tel: 6123263284 Fax: 6122933841 MADE IN CHINA

Bioworld technology, co, Ltd. No 9, weidi road Qixia District Nanjing, 210046, P, R.China. Email: <u>info@biogot.com</u> Tel: 0086-025-68037686 Fax:0086-025-68035151