

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



# TBS Powder

Cat No.: BD5029

### Introduction

This product is Tris Buffered Saline powder; it contains Tris and NaCl. The pH=7.6±0.1 (0.01M), when it is dissolved in purify water. The solution is used for Immunohistochemistry, Immunohistochemistry and Western Blot. It has not been detected for cell culture.

### Reagents

TBS Powder for 2L.

### Application

For washing in IHC, ICC, WB and IF.

### Storage & Shelf life

Store at RT.. Each component is stable for up to 12 Months.

### Procedure

Dissolve it in 2L purify water for use.

### Research Use

For research use only, not for use in diagnostic procedures.

**MADE IN CHINA**

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA. Email: [info@bioworld.com](mailto:info@bioworld.com)

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China. Email: [info@biogot.com](mailto:info@biogot.com)

Tel: 0086-025-68037686

Fax: 0086-025-68035151