

# PRODUCT DATA SHEET

Bioworld Technology, Inc.

# Goat Anti-Mouse IgG (H&L) Highly Cross-Adsorbed - AF488

Catalog: BS22351 Host: Goat

#### **Description:**

Goat Polyclonal Secondary Antibody to Mouse IgG (H&L) Highly Cross-Adsorbed and conjugate to AF488

### **Product:**

1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

#### **Preservative:**

0.02% Proclin/Sodium Azide

#### **Classification:**

Secondary Antibodies

## **Purification&Purity:**

Goat Polyclonal Secondary Antibody to Mouse IgG (H&L) have been cross-adsorbed against IgG from bovine, goat, horse, rabbit and human. Cross-adsorption or pre-adsorption is a purification step to increase specificity of the antibody resulting in less background staining and cross-reactivity. The secondary antibody solution is passed through a column matrix containing immobilized serum proteins from potentially cross-reactive species. Only the nonspecific-binding secondary antibodies are captured in the column, and the highly specific secondaries flow through. Further passages through additional columns result in highly cross-adsorbed preparations of

secondary antibody. The benefits of these extra steps are apparent in multiplexing/multicolor-staining experiments where there is potential cross-reactivity with other primary antibodies or in tissue/cell fluorescent staining experiments where there may be the presence of endogenous immunoglobulins.

## **Applications:**

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

## Storage&Stability:

Shipped at 4  $^{\circ}$ C. Upon delivery aliquot and store at -20  $^{\circ}$ C for one year. Avoid freeze/thaw cycles.

## **Specificity:**

By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins.

## **Expiration Date:**

Six months from date of reconstitution.

#### Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: <u>info@bioworlde.com</u>

Tel: 6123263284 Fax: 6122933841 Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: <u>info@biogot.com</u>
Tel: 0086-025-68037686
Fax: 0086-025-68035151