

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



# WB Blocking Buffer ( phosphorylated antibody specific )

Cat No.: WBA2400-R2

### Introduction

This product is used to block PVDF or NC membranes. The main active can effectively block the protein binding sites on membranes and reduce non-specific bindings of antibodies, thus reducing the background and enhancing the signal-to-noise ratio. The product is suitable for use when detecting phosphorylated antibodies.

### Quantity

100ml/500mL

### Procedure

This product is a ready-to-use liquid and can be directly used for blocking Western Blot membranes. Usually, blocking for 60 minutes at room temperature is sufficient.

### Storage & Shelf life

Store at room temperature. Each component is stable for up to 6 Months.

### NOTE

If the buffer color and smell is changed, you'd better to discard it.

### Research Use

For research use only, not for use in diagnostic procedures.

**MADE IN CHINA**

Bioworld Technology, Inc.  
1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA. Email: [info@bioworld.com](mailto:info@bioworld.com)  
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.  
No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China. Email: [info@biogot.com](mailto:info@biogot.com)  
Tel: 0086-025-68037686 Fax: 0086-025-68035151